

RAPC Chamber

Unpacking Instructions

Revision 1.0

8/19/2004



Bridgekey Corporation
2000 Winton Rd S #5-103
Rochester, NY 14618

(585) 240-6012
www.bridgekey.com

Chamber Unpacking

1. Open the outer shipping carton on top where shipping labels are attached.
2. Remove some of the peanuts and the upper foam bar.
3. Remove the inner shipping carton and place it on a flat work area.
4. Open the inner shipping carton on top as indicated.

Caution

Cut only packing tape!
The top of the chamber is directly below the box
cardboard.

5. Open the inner box flaps and remove the two cardboard end supports.

Caution

The chamber bottom and top are not attached !
When moving the chamber, lift it by holding it under
the bottom plate.

6. Remove chamber from inner carton and place it on a flat work area.
7. Remove outer plastic cover.
8. Remove the black tape holding the top to the chamber.
9. Remove the top.
10. Locate four thumb screws (supplied with the controller and the cables shipment).
11. Install the thumb screws in the lower corners of the chamber to secure the bottom.
12. Remove all black tape from the chamber.
13. Replace top.
14. Chamber is ready for use.

See the RAPC User's Guide for detailed information on the use and maintenance of the chamber, and the installation and use of the software.

Chamber Cleaning

1. Clean the chamber only with water based cleaners, such as Roccal D.
2. Remove the four base thumb screws to separate the bottom from the chamber.
3. Lift the chamber off its bottom plate and set behind the plate.

Caution

Do not use any solvents to clean the chamber !

The chamber acrylic will crack and fracture if it comes into contact with solvent based cleaners or any solvent based substance, such as alcohol.

4. Clean the bottom plate.
5. Set the chamber on its back wall and clean the dividers, if needed.
6. Replace the chamber back onto its bottom plate.
7. Note that the chamber is located on the bottom plate via the four small corner blocks.
8. Replace the four base thumb screws.

For more information on Roccal D, please visit the RAPC web site at www.rapc.us and go to the LINKS page.

Returning The Shipping Boxes

For certain installations, Bridgekey has arranged for return of the shipping boxes and provided prepaid return shipping labels. Please follow the instructions below to repack the boxes and ship them back to Bridgekey.

1. Remove more peanuts from the outer box.
2. Place the empty inner carton back in the outer box.
3. Replace top foam support.
4. Replace all peanuts.
5. Seal the outer box using suitable shipping tape.
6. Apply the return shipping labels.
7. Contact your shipping agent and ship the box back to Bridgekey.
8. Advise Bridgekey of the shipment.

Disposing or Storing The Shipping Boxes

The shipping boxes should be stored if the need for future long distance shipping of the chamber is anticipated. If you do not want to store the boxes or do not have the space for them, please dispose of the shipping materials properly. The peanuts used are biodegradable.

Bridgekey can provide new shipping boxes if the need arises in future.